



TECHNICAL SERVICE BULLETIN # 5

Product: Fabius GS

Date: February 25, 2003

Re: Recall - Phase 2: APL Valve Crater Replacement

Reference Doc: SP00256 – APL Valve Crater Replacement on COSY II

Reason/Rational: Eliminate Possible APL Valve Crater Leaks

Solution/Procedure: Install new APL valve crater assembly via SP00256 on current Fabius GS installed base units. Complete the Return Form attached to SP00256 and return to DMI within 7 days of account completion.

Devices affected:

10001,10002	11372 thru 11377
10004 thru 10448	11379 thru 11385
10450 thru 10631	11388,11389,11396
10633 thru 10645	11398, 11403,11404
10647 thru 10699	11407 thru 11430
10701 thru 10706	11432, 11434 thru 11468
10708 thru 11367	11471 thru 11534
11370	11540 thru 11549
	11552 & 11553

Timing/When: This is the permanent solution to the APL valve crater oxidation. Immediate action is required to prevent breathing system leaks associated to crater valve oxidation. Complete this action ASAP, but no later than 6-months from date of this TSB.

Additional Info: Tools Required:

- 1 - Wrench-Spanner Insert/Extract – P/N 4118165
- 1 - Scalpel – P/N 4118169



Cost:	Valve crater kits come at no charge to the customer. Labor – Warranty
Ordering Info/Parts:	<p><u>US Market</u></p> <p>Valve crater kits are automatically shipped from DMI to Service Representatives (direct) and Distributors according to the attached installed base.</p> <p><u>Non-US Market</u></p> <p>DrägerService worldwide & it's Authorized Service Organizations shall order the appropriate quantity of tools and valve crater update kits necessary to complete this field action on Fabius GS units in their region.</p> <p>Kit; Crater Insert – P/N 4118168 1 - Wrench-Spanner Insert/Extract – P/N 4118165 1 - Scalpel – P/N 4118169</p>
Distribution:	DrägerService & All Authorized Service Organizations worldwide.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047. Fax; 1-215-721-5789 or techsupport@draegermed.com.

